



### **PANDA SCANNER ALL OVER THE WORLD**

⊕ www.panda-scanner.com⋈ info@panda-scanner.com

Chengdu | Ziyang | Suzhou | Ningbo | Hong Kong | Brazil | USA



Create Maximum Value For Every User











# **PARAMETERS**











# FEATURES

## **Product Features**

- Al Technology, much smarter, improve a lot of scanning experience
- 3 Different tips, meet more demands of different applications
- Much smaller, very easy to operate
- (HD) Clear margin line, high accuracy
- (a) Fast scanning speed, easy to learn and use
- (E) True color scan and HD display







## **Application Features**

- High resolution and true image Reproduce the real state of teeth Accurate anatomical features 3D high-definition color, visual communication Easy and exciting treatment experience
- Various range of software Clear and easy Margin Line drawing Pre-op Scan available for Implant and restoration Comparison of two sets of data such as denture and full arch Enhanced scanning option for Veneer
- Case selection and finalized report Auto Preparing the case report to the labs Various options to select in reports such as Bridge or Crown Convenient tools for clinics and laboratories
- Multiple platforms to access the case 3D preview allows to be accessed via the phone Share the cases to patients immediately after scanning Fully prepared case reports to the labs One click link directs to EXOCAD

### **Software Features**





### Undercut Check

Auxiliary confirmation in place, colour gradient indicates undercut, accurate prediction, and improves the communication efficiency between doctors and technicians.

### Margin Line Drawing

Draw accurate and clear margin lines, linking doctors and technicians to reduce rework.

### Pre-preparation Scan

Obtain pre-preparation tooth state, retain bite data, and reshape the original good teeth.

### **Bite Check**

Color gradient indicates distance, accurately judges the restoration distance and assists in improving the effect of tooth preparation.





models.





Full-cycle diagnosis and treatment monitoring, comparative observation of model differences, simultaneous analysis, comprehensive upgrade of diagnosis and research.

Locks the tooth position to prevent data loss and reject excess data interference.

### Al Intelligent Scan

Easily identify soft and hard tissues, eliminate redundant data, and make accurate dental

### **Orthodontic Simulation**

Generate simulation animations to help doctors and patients communicate, share quickly, and promote orders.

### **Disease Course Monitoring**

# **APPLICATIONS**

## Restoration

Accurate and clear shoulder margin brings efficient design, and highdefinition color images help dentists effectively distinguish between gingiva and teeth. One visit, immediate restoration.





## Implant

Fast scanning with a large field of view, easily capture 3mm data of cuff, and accurately scan of metal path pin. No need to do repeat impression and improve patient treatment experience.



## Orthodontics

High accuracy of the full dentition, restore the real state of full arch. Get orthodontic treatment quickly, and save time for more patients.









# **SPECIFICATIONS**

# PRODUCT INFORMATION

Category	PANDA P2
Acquisition Technology	Staring Scan
Camera Number	х З
Scan Field	18x16mm
Scan Depth	0-15mm (Adjustable up to 20mm)
Precision	≤10µm
Accuracy	≤15µm
True Color	Full HD
Anti Fog System	Intelligent Heating
Full Jaw Scanning Time	1-2 min
Handpiece Enclosure	Aviation Aluminum Alloy
Handpiece Dimension	216 x 40 x 36 mm
Handpiece Weight	246g
Tip Types	Stypes (N/M/D)
Number of Tips included	5
Autoclave Cycle for Tips	20-50 times
Calibrator	Automatic
Image Transfer Interface	USB3.0
Cable Lenght (m)	2m
Cart Touchscreen	Optional
Power Supply Type	AC/DC Medical Power Adapter
Supply Voltage (V)	100-240V/50-60Hz
Supply Current (A)	0.7-1.5A
Storage Temperature (°C)	-10°C- 55°C
Operating Temperature (°C)	15°C-30°C

# PC REQUIREMENTS

Category	PAN
CPU	Intel
RAM	16G/3
Hard disk	SSD !
GPU	RTX2
Operating system	Wind
Display	Reso

\*Check the PC Recommendation Catalog on our website for more detail.





### NDA P2

- l i7-12700H/Intel i7-11800H or above
- /32G
- 512G or above
- 2060/RTX3060 or above
- dows10/11 64bit
- olution 1920×1080 and above

# SOFTWARE

Category	PANDA P2
Treatment for	Orthodontics, Implant, Restoration, etc
Image Store	Yes
3D Data Preview	Yes
Output Files	STL, PLY, PTY (encrypted for user)
Cloud Service	Yes
Cloud Service Price	5GB/free
Direct Compatibility with Dental CAD	Exocad*
Main Image Processing Tools	Trim/Edit/Recover/Lock, etc.
Margin Line Drawing	Yes
Occlusal Matching	Yes
Undercut	Yes
Comparison	Yes
AI Technology	Auto Tissue Remove, Overlaping Correct, etc.
License	Free
Software Upgrade	Free
API	Available
Scanning Feedback	Music

\*All names, logos, brands and other trademarks featured or referred to within this brochure are the property of their respective trademark holders.

# AFTERSALES

Category	PANDA P2
Standard Warranty	l year
Extended Warranty	2-3 years available



